

ORDERING

SECTION CONTENTS

- Freedom Product Order Form
(For Manufacturer Direct Orders)
- Component Order Form
- Terms of Sale, Returns & Warranties
- Return Checklist



ORDER FORM - CANADA
WWW.ORTHOACTIVE.COM

TOLL FREE: 1.800.663.1254
FAX: 1.888.558.6565

Order Date: _____ Req. Ship Date: _____ **Ship Via:** Ground 2nd Day Next Day Early AM

Ordered By: _____ PO #: _____ **Ship Notes:** _____

Special Order Notes: _____

BILL TO:

SHIP TO:

Name _____

Name _____

Address _____

Address _____

City _____

City _____

State Zip _____

State Zip _____

Phone () _____

Phone () _____

Note: If Warranty or 60 Day Return Re-Order, please call Customer Service prior to ordering.

Patient Data

Name / ID#: _____ Age: _____ Gender: M F

Foot Size (cm): _____ Weight (lbs): _____ Height: _____ Ft _____ in _____

Amputation: Left Right Bilateral Level: Transtibial - BK Transfemoral - AK Hip Disartic - HD Knee Disartic - KD Symes

Module Needed: Left Right Both Patient Notes: _____

Skin Tone: Beige Brown _____

_____ Low Activity (K2 or K3) Typical daily living, walking, golf, bicycling

_____ Moderate Activity (K3) Repetitive activities, recreational sports, jogging, dancing

_____ High Activity (K3 to K4) Rigorous activities, competitive sports, running, power lifting, snowboarding

Clearance Measurement (inches): _____ (From most distal aspect of socket to the ground
Include space for a liner, shuttle lock, etc.)

RETURNS, ADJUSTMENTS, CREDITS


We are committed to the complete satisfaction of the prosthetist and amputee. We deliver most products with a 60 day "Satisfaction Guarantee." If you are dissatisfied with a Freedom standard production product for any reason, you may return it for a full refund within 60 days of the original invoice date. Modular components and liners can be returned within 60 days of purchase only if their packaging is unopened. All returns require a Return Authorization Number (RA#), which is obtained by calling Freedom Innovations Customer Service or the Distributor through which the product was originally ordered. No returns will be accepted without an RA #.

A restocking fee of 20% may be applied to orders returned after the 60 day Satisfaction Guarantee.

Internal Use CS Rep: _____ Date: _____ Order #: _____
MF Rep: _____ Notes: _____

PRODUCT SELECTION- Page 2


(Please complete appropriate data fields for your selection.)




Renegade



Renegade LP




Runway



Silhouette



Silhouette LP



Nitro



Slalom Ski

- Renegade® (RS1)**
- Renegade® LP (RS6)**

Sizes (22-30 cm) _____

Stiffness Category (1-9) _____

Foot Shell Height Mid-profile (RS1 Default) High-profile (RS6 Default) No Shell

Optional Pre-Cut Shin Foam

Special Order Size 31 cm (7-10 day lead time)

- Renegade® SX (RS1-S1) Single-Axis Rotation**
- Renegade® LP-SX (RS6-S1)**

Sizes (22-30 cm) _____

Stiffness Category (1-7) _____

Foot Shell Height SX Mid-profile (RS1-S1 Default) High-profile (RS6-S1 Default) No Shell

Special Order Size 31 cm (7-10 day lead time)

- Renegade® MX (RS1-M1) Multi-Axial Rotation**
- Renegade® LP-MX (RS6-M1)**

Sizes (22-30 cm) _____

Stiffness Category (1-7) _____

Foot Shell Height MX Mid-profile (RS1-M1 Default) High-profile (RS6-M1 Default) No Shell

Special Order Size 31 cm (7-10 day lead time)

- Runway® (RS2) (Push button adjustable)**
- Runway® HX (RS2-K1) (Hex key adjustable)**

Sizes (22-28 cm) _____

Stiffness Category (1-6) _____

Foot Shell Style Standard Sandal Toe No Shell

- Silhouette™ (RS4) w/34mm female tube clamp**
- Silhouette™ LP (RS5) w/male pyramid adapter**

Sizes (22-30 cm) _____

Stiffness Category (1-9) _____

Foot Shell Height Mid-profile (Default) High-profile No Shell

Special Order Size 31 cm (7-10 day lead time)

30 mm knee pyramid adapter

Versa Clamp

CUSTOM PRODUCTS

Nitro Running Foot (7-10 day lead time)

User Weight (lbs.) _____

Tall Profile (9.5") (FS6) Short Profile (7.5") (FS9)

Slalom (SKI) (7-10 day lead time)

Length of Sound Side Ski Boot in cm _____

OR Length of Ski Binding in cm _____

User Weight (lbs.) _____

Freestyle Swim Foot (SW2) (24 hour lead time)

Sizes (22-30 cm) _____

Stiffness Category (1-6) _____

Foot Shell Height Mid-profile No Shell

Special Order Size 31 cm (7-10 day lead time)

PRODUCT INFORMATION- Page 4

| Product Clearance & Weight Chart | | | Avg. Weights | |
|--|----------|-----------|--------------|-------|
| Product | Size | Clearance | Foot | Shell |
| Sierra (FS1) | 22-25 cm | 5.2 in. | 343g | 127g |
| | 26-28 cm | 6.0 in. | 404g | 159g |
| | 29-31 cm | 6.7 in. | 493g | 219g |
| Highlander (FS3) | 22-25 cm | 3.7 in | 310g | 127g |
| | 26-28 cm | 4.2 in | 371g | 159g |
| | 29-31 cm | 4.7 in | 461g | 219g |
| Pacifica (FS2) | 24-25 cm | 2.1 in | 259g | 127g |
| | 26-28 cm | 2.2 in | 310g | 159g |
| | 29-31 cm | 2.4 in | 382g | 219g |
| Pacifica LP (FS4) | 22-25 cm | 1.8 in | 305g | 127g |
| | 26-28 cm | 2.0 in | 370g | 159g |
| | 29-31 cm | 2.2 in | 435g | 219g |
| Renegade (RS1) | 22-25 cm | 5.8 in | 358g | 127g |
| | 26-28 cm | 6.6 in | 437g | 159g |
| | 29-31 cm | 7.3 in | 549g | 219g |
| Renegade LP (RS6) | 22-25 cm | 4.4 in | 337g | 127g |
| | 26-28 cm | 4.8 in | 403g | 159g |
| | 29-31 cm | 5.2 in | 496g | 219g |
| Renegade SX (RS1-S1) | 22-25 cm | 8.4 in | 765g | 127g |
| | 26-28 cm | 8.5 in | 836g | 159g |
| | 29-30 cm | 8.6 in | 907g | 219g |
| Renegade LP-SX (RS6-S1) | 22-25 cm | 6.8 in | 632g | 127g |
| | 26-28 cm | 6.9 in | 725g | 159g |
| | 29-30 cm | 7 in | 794g | 219g |
| Renegade MX (RS1-M1) | 22-25 cm | 7.2 in | 672g | 127g |
| | 26-28 cm | 7.25 in | 742g | 159g |
| | 29-30 cm | 7.4 in | 794g | 219g |
| Renegade LP-MX (RS6-M1) | 22-25 cm | 5.6 in | 553g | 127g |
| | 26-28 cm | 5.65 in | 646g | 159g |
| | 29-30 cm | 5.7 in | 717g | 219g |
| Runway (RS2) w/Mid-Profile Shell | 22-25 cm | 3.0 in | 367g | 158g |
| | 26-28 cm | 3.25 in | 449g | 210g |
| Runway/HX (RS2-K1) w/Sandal Toe Shell | 22-25 cm | 3.0 in | 372g | 131g |
| | 26-28 cm | 3.25 in | 414g | 414g |
| Silhouette (RS4) | 22-25 cm | 7.2 in | 423g | 127g |
| | 26-28 cm | 7.7 in | 453g | 159g |
| | 29-31 cm | 8.1 in | 483g | 219g |
| Silhouette LP (RS5) | 22-25 cm | 6.8 in | 357g | 127g |
| | 26-28 cm | 7.3 in | 387g | 159g |
| | 29-31 cm | 7.7 in | 417g | 219g |
| LP Symes (LP2) | 22-25 cm | 1.7 in | 227g | 127g |
| | 26-28 cm | 1.9 in | 295g | 159g |
| | 29-31 cm | 2.1 in | 363g | 219g |
| Senator (VS1) | 22-25 cm | 4.5 in | 362g | 127g |
| | 26-28 cm | 5.0 in | 390g | 159g |
| | 29-30 cm | 5.4 in | 418g | 219g |
| Defender (FS3) | 16-18 cm | 2.4 in | 92g | 80g |
| | 19-21 cm | 3.0 in | 107g | 121g |
| Freestyle (LP2-W2) | 22-25 cm | 3.99 in | 485g | 127g |
| | 26-28 cm | 4.07 in | 542g | 159g |
| | 29-30 cm | 4.55 in | 635g | 219g |

Determine Your Impact Level

| Low Hobby Oriented Athlete | Moderate Recreational Athlete | High Extreme Athlete |
|-------------------------------|----------------------------------|-------------------------|
| Walking | Jogging | Baseball |
| Gardening | Tennis | Running |
| Bowling | Volleyball | Power Lifting |
| Shopping | Hiking | Soccer |
| Golfing | Dancing | Football |
| Biking | Skiing | Skateboarding |
| Fishing | Martial Arts | Snowboarding |
| Camping | Boxing | Competitive Sports |
| Surfing | Baseball | |
| | Skating | |

Stiffness Category Selection Chart

| Weight (lbs) | Category | | | Weight (kg) |
|--------------|----------|----------|------|-------------|
| | Low | Moderate | High | |
| 100-115 | 1 | 1 | 2 | 44-52 |
| 116-130 | 1 | 2 | 3 | 53-59 |
| 131-150 | 2 | 3 | 4 | 60-68 |
| 151-170 | 3 | 4 | 5 | 69-77 |
| 171-195 | 4 | 5 | 6 | 78-88 |
| 196-220 | 5 | 6 | 7 | 89-100 |
| 221-255 | 6 | 7 | 8 | 101-116 |
| 256-285 | 7 | 8 | 9 | 117-130 |
| 286-325 | 8 | 9 | | 131-147 |
| 326-365 | 9 | | | 148-166 |

Senator Stiffness Category Chart

| Weight (lbs) | Category | | Weight (kg) |
|--------------|----------|--------|-------------|
| | Low | Medium | |
| 100-120 | 1 | 2 | 45-55 |
| 121-145 | 2 | 3 | 56-66 |
| 146-175 | 3 | 4 | 67-80 |
| 176-210 | 4 | 5 | 81-96 |
| 211-250 | 5 | 6 | 97-114 |
| 251-300 | 6 | | 115-136 |

Defender Stiffness Category Chart

| Weight (lbs) | Category | Weight (kg) |
|--------------|----------|-------------|
| 20-30 | 1P | 9-13 |
| 31-40 | 2P | 14-18 |
| 41-50 | 3P | 19-22 |
| 51-60 | 4P | 23-27 |
| 61-70 | 5P | 28-31 |
| 71-80 | 6P | 32-36 |
| 81-100 | 7 Adult | 37-45 |

























Distributed in Canada
by Ortho Active
Toll Free: 1.800.663.1254
Fax: 1.888.558.6565

COMPONENT ORDER FORM

| | | |
|---------|-------------|--------------|
| Sale #: | Ordered By: | Customer PO# |
|---------|-------------|--------------|

| | |
|-----------------|-----------------------|
| BILL TO: | SHIP TO: |
| Name/Company | Name/Company |
| Address | Address |
| City, ZIP | City, ZIP |
| () () | Telephone/FAX () () |

| Qty | Component | | Qty | Component | | Qty | Component |
|-----|--|--|-----|---|--|-----|---|
| | AC-125 30mm Female Tube Clamp  | | | AC-130 4-hole Male Pyramid  | | | AC-160 4-prong Male Pyramid Base  |
| | AC-126 30mm Male Tube Clamp  | | | AC-132 4-hole Female Pyramid  | | | AC-161 4-hole Lamination Base  |
| | AC-127 30mm x 100mm Titanium Pylon w/ Female Pyramid  | | | AC-137 4-hole Pyramid Receiver  | | | AC-170 32mm Dual Female Adaptor  |
| | AC-128 30mm x 400mm Aluminum Pylon  | | | AC-138 3-prong Pyramid Receiver  | | | AC-171 45mm Dual Female Adaptor  |
| | AC-145 30 x 3mm x 15" Carbon Pylon  | | | AC-139 Rotating Male Pyramid  | | | AC-172 60mm Dual female Adaptor  |
| | AC-146 30 x 3mm x 9" Carbon Pylon  | | | AC-140 Rotating Female Pyramid  | | | AC-173 75mm Dual Female Adaptor  |
| | AC-148 30 x 4mm x 15" Carbon Pylon  | | | AC-135 4-hole Receiver with Male Pyramid  | | | AC-101-18 18 inch Spectra Sock |
| | AC-149 30 x 4mm x 9" Carbon Pylon  | | | AC-136 4-hole Receiver with Female Pyramid  | | | AC-101-12 12 inch Spectra Sock |

Requested Arrival _____ Ship Method [] Next Day [] 2nd Day or [] _____

Notes: